

Flexible

Patented and Patents Pending
See watcomfg.com for details.

Innovator® Flex Series

Flexible Tubing • Built-In Testing & Snap-On Overflow Plate • Coarse Thread: 11.5 tpi

WatcoFlex® Direct Drain Only

Complete Bath Wastes
Rough-In Kits
Stopper/Faceplate Kits
PVC and ABS Fittings



WITH OPTIONAL W-OVERFLOW COVER

Fits Kohler's Tea-for-Two tub.



FLEXW ___ - ___ - CP
BN-BZ-CB-NP-PB.....

Each WatcoFlex® bath waste is custom-fit to each bathtub model.

PLEASE PROVIDE THE TUB MANUFACTURER AND MODEL NUMBER WHEN ORDERING.

FLEX___: The length portion (___) of the WatcoFlex® model number will be determined by the tub's dimensions.

Choose either a complete bath waste OR a rough-in kit and stopper/faceplate kit.

Example Product #: $\overset{u}{\text{FLEX}}\overset{v}{230} - \overset{w}{\text{PF}} - \overset{x}{\text{PVC}} - \overset{y}{\text{CP}} - \overset{z}{***}$

Series Number Tub Closure Material Finish Option

12 PER CARTON

LT

LIFT AND TURN

**COMPLETE BATH WASTES
LENGTHS 14.5" TO 17" LONG**

FLEX ___ - LT - ___ - CP
FLEX ___ - PP - ___ - CP
FLEX ___ - PF - ___ - CP

PF: NOT in BO-BS-WH

FLEX ___ - FA - ___ - CP
FA: BN-BZ-CB-CP-NP-PB

Other Finishes:

AP	BS	PB
BB	BZ	WH
BN	CB	WI
BO	NP	

**RI (ROUGH-IN)
FITTINGS - TUBING - STRAINER BODY
LENGTHS 14.5" TO 17" LONG**

FLEX ___ - LTRI - ___ - CP
FLEX ___ - PPRI - ___ - CP ...
FLEX ___ - PFRI - ___ - CP ...

PF: NOT in BO-BS-WH

FLEX ___ - FARI - ___ - CP
FA: BN-BZ-CB-CP-NP-PB

Other Finishes:

AP	BS	PB
BB	BZ	WH
BN	CB	WI
BO	NP	

**900-SERIES SF (Flex ___)
STOPPER/FACEPLATE
NOT A TRIM KIT - NO STRAINER BODY**

900-LTSF-0-CP
900-PPSF-0-CP
900-PFSF-0-CP

PF: NOT in BO-BS-WH

900-FASF-0-CP
FA: BN-BZ-CB-CP-NP-PB

Other Finishes:

AP	BS	PB
BB	BZ	WH
BN	CB	WI
BO	NP	

PP

PUSH PULL®

**COMPLETE BATH WASTES
LENGTHS 17.5" AND LONGER**

FLEX ___ - LT - ___ - CP
FLEX ___ - PP - ___ - CP
FLEX ___ - PF - ___ - CP

PF: NOT in BO-BS-WH

FLEX ___ - FA - ___ - CP
FA: BN-BZ-CB-CP-NP-PB

Other Finishes:

AP	BS	PB
BB	BZ	WH
BN	CB	WI
BO	NP	

**RI (ROUGH-IN)
FITTINGS - TUBING - STRAINER BODY
LENGTHS 17.5" AND LONGER**

FLEX ___ - LTRI - ___ - CP
FLEX ___ - PPRI - ___ - CP ...
FLEX ___ - PFRI - ___ - CP ...

PF: NOT in BO-BS-WH

FLEX ___ - FARI - ___ - CP
FA: BN-BZ-CB-CP-NP-PB

Other Finishes:

AP	BS	PB
BB	BZ	WH
BN	CB	WI
BO	NP	

**900-SERIES SF (Flex ___)
STOPPER/FACEPLATE
NOT A TRIM KIT - NO STRAINER BODY**

900-LTSF-0-CP
900-PPSF-0-CP
900-PFSF-0-CP

PF: NOT in BO-BS-WH

900-FASF-0-CP
FA: BN-BZ-CB-CP-NP-PB

Other Finishes:

AP	BS	PB
BB	BZ	WH
BN	CB	WI
BO	NP	

PF

PRESFLO®
ADA

FA

FOOT ACTUATED
ADA

QuickTrim® Kits



Innovator® QuickTrim Kits
See page 52

Install chrome on the rough-in and easily upgrade to a special finish without removing the strainer body! See pages 48-49.

- No reduction in flow rate.
- Produces very little additional height or inside-diameter reduction.
- Uses an o-ring for a tight, snug fit (one for brass strainer bodies and one for zinc strainer bodies).
- Eliminates destroying the integrity of your test and the possibility of cross-threading or damaging the tub when replacing the strainer body.